



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6707



Rapid Dissolve Tablet Mold (RDT), Profiller

96 Cavities, 0.940 mL/cavity

Description: This device can be used to make 96 rapid dissolve tablets simultaneously and can also be used for commercially available bases.

This mold is autoclave safe.

Mold Cleaning and Maintenance:

1. Rinse the cavity plate under running lukewarm tap water. If necessary, soak the mold in warm water-soap solution for few minutes to facilitate cleaning.
2. Wash the molds using warm water, liquid soap and a soft brush or sponge. Do not use acidic or caustic cleaning solutions.
3. Rinse and drain off cavity mold and tamper plate.
4. Use 50% alcohol for final wipe down.
5. Let dry. May be dried in a hot air oven.

ATTENTION : Avoid use of metal utensils (spoon, spatulas, scrapper, brush, etc.) and cleaning aids with abrasive surface at any time due to potential risk of damaging coating of the molds.

ProFiller Mold can be cleaned in dishwasher on rinse cycle mode only without the use of dishwasher detergents and with all the above precautions.

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

TECHNICAL SPECIFICATIONS:

Overall Dimensions, W × D × H (mm / in)	235 × 330 × 205 / 9.3 × 13.0 × 0.8	
Weight, Tamper and Cavity Plate (kg / lbs)	3.35 / 7.38	
Individual Cavity, Ø × H (mm / in)	13.4 × 7.5 / 0.53 × 0.30	
Hole volume (mL)	0.940	
Material	Tamper Plate	Hard anodized aluminum
	Cavity Plate	Teflon coated Hard anodized aluminum

USA | CANADA | AUSTRALIA PACIFIC | EUROPE
WWW.MEDISCA.COM