



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.:5731



Hood Accessories, SafeGuardPlus 1200 Trolley

DESCRIPTION: This base stand provides mobility to your hood setup and is equipped with a lower storage shelf. Locking casters fix the hood in place.

APPLICABLE PRODUCTS:

 5727: Hood, SafeGuardPlus 1200, Double Hepa Filter, H985mm x W1220mm x D620mm

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM