




New & Improved Medi-RDT™ Mold



You Asked. We Listened.

Now available in highly durable hospital-grade polymer, PPSU, the Medisca Medi-RDT Mold is designed for repeated and prolonged use. The newly designed tamper with two handles on each side is also advantageous, allowing an even tamp without adding extra stress to the tamper and mold.



Features & Benefits

- Teflon-Free with no resurfacing required.
- Made of hospital-grade polymer, PPSU.
- 96 cavities with a volume of 0.93 mL, yielding an average tablet weight of 750 mg using the Medisca Medi-RDT™ Base.
- Designed for repeated and prolonged baking, up to 120 °C (248 °F).
- Lateral containment system that helps avoid powder loss.
- Lightweight, easy handling, and cleaning. (Mold is dishwasher & autoclave safe; tamper is not)




**Scan QR Code to
BUY NOW**

Product No.: 7314


Medi-RDT™ Kit

Packaged into a Kit Format for Your Convenience



Scan QR Code to
BUY NOW

Product No.: 7218



Formulas

Formula No.	APIs and Dosages	Base
F.005.468v5	Progesterone 100 mg Oral Rapid-Dissolve Tablets	Medi-RDT Base
F.006.495v2	Biotin 5 mg Oral Rapid-Dissolve Tablets	Medi-RDT Base
F.005.470v4	Testosterone 20 mg Oral Rapid-Dissolve Tablets	Medi-RDT Base
F.005.472v4	Melatonin 5 mg Oral Rapid-Dissolve Tablets	Medi-RDT Base

For More Information

USA 1-800-932-1039 medisca.com

