

PRODUCT NO.: 5403



Tuttnauer T-Edge 11 Autoclave, Class B (110V)

Description: Fully automatic table top autoclave that uses steam as a sterilizing agent.

Features:

- Touchscreen used for monitoring and control purposes
- Built -in USB port to enable the operation of an external optional barcode printer
- Built-in memory to record up to 999 sterilization cycles
- Electrical door locking pin that blocks door opening during operation
- Safety thermostat, to prevent over-heating of the chamber heating elements

Technical Specifications:

External Dimensions (Width x Height x Depth)	19.7" x 19.7" x 22.8"
Chamber Dimensions (Diameter x Depth x Volume)	11" x 18" x 27.2 L
Max. Allowable Working Pressure	~40.6 PSI
Tray Dimensions	16.6" x 8.1"x 0.8"
No. of Trays	5
Net Weight	124 lbs / 56 kg
Shipping Weight	152lbs / 69 kg
Total Power	1500W
Voltage	110V
Amperage	12.5A
Frequency	60Hz

USA | Canada | Australia | International

medisca.com