



Luer Lock to Luer Lock, Female-Female Connector, White (Sterile)

Description: Sterile connector facilitates the transfer (or mixing of fluids) between syringes and the connecting of various IV sets.

Features:

- Individually packed
- Ethylene Oxide (EO) sterilized
- DEHP, BPA & latex free
- CE marked

Technical Specifications:

Material	Polypropylene
Height (mm / in)	29.2 / 1.2
Diameter (mm / in)	12.7 / 0.5

USA | Canada | Australia | International

medisca.com