



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 5516**



## 6' Vertical Laminar Flow Cabinet (Purair LF Series)

**Description:** This 6' Vertical Laminar Flow (VLF) hood is ideally suited for use with non-hazardous contaminants.

### FEATURES:

- ✔ UPLA main filter, 99.999% efficient at particle size 0.12 µm
- ✔ Pre-Filter made of disposable polyester fibers with 85% arrestance
- ✔ Air cleanliness meets and exceeds ISO Class 4
- ✔ White epoxy coated steel frame exterior finish
- ✔ Low noise and vibration levels
- ✔ Stainless steel work surface
- ✔ Frame and work surface is durable and chemically resistant
- ✔ Compliant with criteria set forth by USP 797 for sterile preparation in pharmaceutical compounding of nonhazardous agents




### TECHNICAL SPECIFICATIONS:

External Dimensions (W x D x H)	76.4" x 29.25" x 47.25" / 1941 x 743 x 1200 mm
Shipping Dimensions (W x D x H)	91" x 44" x 54" / 2311 x 1118 x 1372 mm
Nominal Width	72" / 1829 mm
Internal Height	28.25" / 718 mm
Internal Depth	28.25" / 718 mm
Net Weight	478 lbs / 217 kg
Ship Weight	595 lbs / 270 kg
Airflow	0.45 m/s (90 fpm)
Electrical	120V, 60Hz or 230V, 50Hz voltages available

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)

## Related Product / Additional Accessories:

-  Hood Accessories, 6' Base Stand w/Lockable Casters (Purair LF Series) – 5517
-  Hood Accessories, 6' IV Bar (Purair LF Series) – 5515
-  Hood Accessories, 6' UV Light (Purair LF Series) – 5514