








PRODUCT NO.: 5702



Syringe, Oral, Clear, 1mL, w/ Tip Cap
(White Plunger)

Description: These single use non-sterile oral syringes have a two-piece design that easily and smoothly administers medication.

FEATURES:

-  Unique tip design reduces residual content
-  Includes self-righting tip caps
-  Class 1 medical device
-  CE marked
-  Latex, BPA and PVC free

TECHNICAL SPECIFICATIONS:

Material	
Barrel	Polypropylene (PP)
Plunger	High density polyethylene (HDPE)
Dimensions	
Length	3-3/4" - 6-1/8" when extended to 1.0 mL
Outside diameter	1/4"
Tip dead space volume	0.045 mL

USA | CANADA | AUSTRALIA PACIFIC | EUROPE
WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.