



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING



PRODUCT NO.: 5810



Filler, Edhard Spout, Blunt End

DESCRIPTION: The spout attaches to the Edhard Extension Assembly and a spout adapter. It is used to fill small diameter tube shaped containers. The length of the Spout reaches to the bottom of the tube.

APPLICABLE PRODUCTS:

-  5811: Filler, Edhard Spout, Adapter **(REQUIRED)**
-  5813: Filler, Edhard Extension Assembly **(REQUIRED)**

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM