




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5811



Filler, Edhard Spout, Adapter

DESCRIPTION: The Adaptor attaches to the Edhard Extension Assembly and a spout. It is used to fill small diameter tube shaped containers. The length of the Spout reaches to the bottom of the tube.

APPLICABLE PRODUCTS:

-  5813: Filler, Edhard Extension Assembly (**REQUIRED**)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM