




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5816







Filler, Edhard Pump Housing Assembly

DESCRIPTION: The pump housing assembly allows to dispense viscous and delicate products, containing particles up to 10mm diameter.

TECHNICAL SPECIFICATIONS:

Pump Configuration	Lobe Type
Impeller	5 lobe type, HDPE

APPLICABLE PRODUCTS:

-  5814: Filler, Edhard Hopper, 4 Qt., Stainless steel **(REQUIRED)**
-  5815: Filler, Edhard Hopper, 6 Qt., Tritan **(REQUIRED)**
-  5817: Filler, Edhard Power Base, 220/240 V, 50/60 Hz **(REQUIRED)**
-  5818: Filler, Edhard Power Base, 120V, 50/60 Hz **(REQUIRED)**

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM



MEDISCA®
YOUR TRUSTED PARTNER IN COMPOUNDING

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.