PRODUCT NO.: 5830





Filter, Syringe, D65R, PTFE Membrane, 65mm, Dual Layer, 1.0µm/0.2 µm (Sterile) Hydrophobic, Female Luer-Lock Inlet, Male Luer-Lock Outlet

DESCRIPTION: These filters are designed for applications where sterility of the affluent is mandatory. The product contains a pharmaceutical grade PTFE membrane (Sterilizing Grade) and the filters are built in a 10,000-class clean room, flushed, dried and 100% integrity tested prior to packaging. The product is factory sterilized by Ethylene Oxide gas (ETO) ensuring sterility.

TECHNICAL SPECIFICATIONS:

Material of Construction	Filter Membrane: PTFE
	Supports/Capsule Body: Polypropylene
Filter Diameter (mm/in)	65
Pore Size (µm)	1.0/0.2
Process Volume (mL)	460-1880
Inlet Connection	Female Luer-Lock
Outlet Connection	Male Luer-Lock
Min. Bubble Point (25°C)	≥ 3.5 bar/50 psi
Max. Operating Temperature	80°C/176°F
Max. Operating Pressure	2.4 bar/35.3 Psi
UPS Bacteria Endotoxins	<25 EU/ml
Sterilization	Ethylene Oxide Gas (ETO)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.