



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5869



Hood Accessories, Exhaust Collar (Purair Advanced)

DESCRIPTION: The exhaust collar allows the enclosure to be thimble-ducted to the outside for additional protection.

TECHNICAL SPECIFICATIONS:

- ☑ Collar Diameter: 8"

APPLICABLE PRODUCTS:

- ☑ 5970: 3' Ductless Fume Hood (Purair Advanced)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM