



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5899



Water Purification System, WaterPro RO, Dispensing Gun

DESCRIPTION: Dispensing Gun provides additional means to dispense water from the WaterPro RO Station. Trigger permits operation with one hand and provides precise drop-by-drop control or locks to provide a continuous stream of water. Includes gun attached to a 6' hose, epoxy-coated steel gun rest and mounting hardware.

TECHNICAL SPECIFICATIONS:

Weight (lbs/kg)	5/2.3
-----------------	-------

APPLICABLE PRODUCTS:

- 5897: Water Purification System, WaterPro RO, Wall Mount
- 5898: Water Purification System, WaterPro RO, Membrane
- 5900: Water Purification System, WaterPro RO, Carbon Prefilter Kit

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM