



PRODUCT NO.: 5941



Hood Accessories, HEPA Filter (Air Science Fume Box)

DESCRIPTION: This filter is for fine powders and particulates. Use alongside the pre-filter to remove contaminants from the air.

TECHNICAL SPECIFICATIONS:

Nominal Depth (in)	9
Nominal Height (in)	1-1/2
Nominal Width (in)	19
Filter Efficiency	99.97% (Particle Size 0.3 Microns)
Shipping Weight (lbs)	25

APPLICABLE PRODUCTS:

-  6942: Hood, 2' Ductless Fume Box, Low Profile (Air Science)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM