




**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING





**PRODUCT NO.: 5948**



## Forceps, Polystyrene w/ Blunt tip (Sterile)

**DESCRIPTION:** These 4.9" long plastic tweezers/forceps feature serrated inside tips and ribbed finger supports.

### FEATURES:

-  Made of polystyrene
-  Serrated blunt tip
-  Disposable
-  Individually peel-packed

### TECHNICAL SPECIFICATIONS:

Length (in/cm)	4.9/124.46
Width (in/cm)	0.9/22.86
Tip opening (in/cm)	0.6/15.24

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)