



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 5958**



## **Hood Accessories, Cart/Stand, for Powdersafe Type B (Airclean), 32"**

**DESCRIPTION:** Metal cart with casters constructed of 14-gauge epoxy coated steel. Carts feature a two-position movable base bar for wheelchair or lab chair compatibility.

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)