



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5988



Hood Accessories, Primary HEPA Filter (Purair Advanced)

DESCRIPTION: HEPA filter for Purair Advance Models. These self-contained filters are designed to physically capture particles larger than 0.3 microns with >99.999% typical efficiency.

TECHNICAL SPECIFICATIONS:

External Dimensions, W x D x H (in)	14 x 20 x 4
Filter Type	HEPA
Material	Plastic, Carbon

APPLICABLE PRODUCTS:

-  5989-01: 6' Ductless Fume Hood (Purair Advanced)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM