



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5994



Hood Accessories HEPA Filter, 42" (Airclean)

4ft, For AirClean Type B PowderSafe Balance Enclosure

DESCRIPTION: Replacement HEPA filter to be used for the 4 ft PowderSafe Balance Enclosure. The filter has a 99.997% efficiency rating for collecting 0.3 micron particles.

APPLICABLE PRODUCTS:

 7563-01: Hood, 4' PowderSafe, Type B Balance Enclosure, W/AirFlow M

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM