



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5995



Cylinder, Graduated, Polypropylene, 1000mL

DESCRIPTION: Polypropylene Graduated Cylinders 1000mL cylinders feature a single scale with raised numeral graduations in 10mL increments. Non-wetting polypropylene surface eliminates concave menisci.

FEATURES:

- ✓ Translucent cylinder
- ✓ Molded-in large graduations
- ✓ Chemical resistant

TECHNICAL SPECIFICATIONS:

Capacity (mL)	1000
Graduations (mL)	10
Subdivison (mL)	10
Base Type	Circular, ribbed
Autoclavable	Yes
Material	Polypropylene (PP)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM