



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5996



Cylinder, Graduated, Polymethylpentene, 100mL

Class A

DESCRIPTION: This Class A certified graduated cylinder meets DIN 12681, ISO 6706 standards for Volumetric Glassware.

TECHNICAL SPECIFICATIONS:

Diameter (in/cm)	1.29/3.3
Height (in/mm)	9.49/250
Capacity (mL)	100
Graduations (mL)	1
Subdivison (mL)	1
Tolerance (mL)	± 0.5
Autoclavable	Yes
Material	Polymethylpentene (PMP)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM