



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5997



Cylinder, Graduated, Polymethylpentene, 50mL

Class A

DESCRIPTION: This Class A certified graduated cylinders meet DIN 12681, ISO 6706 Standards for Volumetric Glassware.

TECHNICAL SPECIFICATIONS:

Diameter (in/cm)	1.06/2.7
Height (in/mm)	7.87/200
Capacity (mL)	50
Graduations (mL)	1
Tolerance (mL)	± 0.5
Material	Polymethylpentene (PMP)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM