



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 5998**



## Cylinder, Graduated, Polymethylpentene, 10mL

**DESCRIPTION:** This Class A certified graduated cylinder meets DIN 12681, ISO 6706 standards for Volumetric Glassware.

### TECHNICAL SPECIFICATIONS:

Capacity (mL)	10
Tolerance (mL)	0.10
Height (in/mm)	5.70/145
Autoclavable	Yes
Subdivision (mm)	0.20
Graduations (mL)	0.2
Max Teperature (°C)	150
Material	Polymethylpentene (PMP)

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)