



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6002



Beaker, Polypropylene, 5000 mL

DESCRIPTION: This reusable polypropylene low-form beaker meets DIN, ASTM Class A and B certifications. This beaker is used to hold and sometimes measure liquids.

TECHNICAL SPECIFICATIONS:

Diameter (in/cm)	7.12/18.1
Height (in/cm)	9.7/24.7
Capacity (mL)	5000
Graduated	Yes
Increments (mL)	500
Color	Translucent
Style	Griffin
Scale	Blue graduated scale
Autoclavable	Yes
Disposable	No
Temperature Max. Operating (°C/°F)	20/68
Type	Beaker
Material	Polypropylene

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM