



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6029



4' Hydraulic Lift Stand w/Fixed Feet, 32-42" Working Height

DESCRIPTION: Electric Hydraulic Lift Base Stands offer the flexibility to be adjusted at any time. Height is infinitely adjustable from 32" to 42". They allow you to position the stand and the enclosure is supports to accommodate persons of different heights and wheelchair-bound users. This is applicable to XPert balance enclosures and Xpert/RXPert filtered balance systems.

TECHNICAL SPECIFICATIONS:

Dimensions, W x D (in/cm)	49.2 x 34.3/125 x 87.1
Handle/Switch location	Right Side
Hydraulic lift	Electric
Maximum Stand Height (in/cm)	42/106.7
Minimum Stand Height (in/cm)	32/81.3
Nominal Width	4'
Shipping Weight (kg/lbs)	79/175
Electrical (V/Hz/A)	115/60/2.4

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM