



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING




**PRODUCT NO.: 6032**



## **Weigh Boat, Anti-Static, 1.6 x 1.6 x 0.3", Small, White**

**DESCRIPTION:** Anti-static disposable polystyrene weight boats are flexible for handling liquids and solids. Square with stable rounded bottom, smooth surface and contoured sides.

### **FEATURES:**

-  Easy pour spout
-  Boats are shallow with a wide, flat bottom to resist tipping
-  Rounded corners simplify transfer

### **TECHNICAL SPECIFICATIONS:**

Capacity (mL)	7
Dimension (mm/in)	41 x 41 x 8/1.6 1.6 x 0.3
Color	White
Material	Polystyrene (PS), Anti-Static