**PRODUCT NO.: 6038** 





## Weigh Canoe, Anti-Static, 7.1 x 4.3 x 1", XLarge, White

**DESCRIPTION:** This weigh canoe protects the pans of weighing balances from damage and provides a convenient way to transfer weighed materials to another container.

## **FEATURES:**

- **2** Pour-spout design easily dispenses samples without flexing the dish
- Z Flat-bottom dishes are from 12-mil thick polystyrene for better strength, stability and handling

## **TECHNICAL SPECIFICATIONS:**

Capacity (mL)	240
Dimension (mm/in)	181 x 110 x 26/7.1 x 4.3 x 1
Color	White
Material	Polystyrene (PS), Anti-Static