




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING



PRODUCT NO.: 6040



Weigh Boat, Hexagonal, Anti-Static, OD:3" Dep:3/4", Medium, White

DESCRIPTION: This hexagonal shape weigh boat provides accurate pouring with minimal sample loss. Wide with flat bottom to resist tipping.

FEATURES:

-  Bend flexible boats at point to form pour spouts
-  Flat bottoms and sloping sides provide accurate transfer

TECHNICAL SPECIFICATIONS:

Capacity (mL)	50
Dimension (in)	OD: 3, Dep: 3/4
Color	White
Material	Polystyrene (PS), Anti-Static