



DISPOSABLE WEIGH BOATS & CANOES



Antistatic polystyrene ensures minimal product loss



Also available in black for efficient color contrasting



Effective cost solution

These flexible, antistatic, polystyrene weigh boats and canoes feature a smooth surface that enables accurate pour-outs with minimal product loss. This is a fast, safe and reliable way to measure and transfer solid or liquid products. MEDISCA weigh boats and canoes are ideal for use with balances, microwaves and moisture analyzers. They are available in transparent light white and black for accurate results.

FEATURES & BENEFITS

Biologically inert and resistant to diluted and weak acids, aqueous solutions, alcohols and bases

Flexible material allows slight folding for easy transfer of product

Temperature resistant from -10°C to 70°C (14°F to 158°F)

Smooth surface with rounded corners simplifies pouring

Shallow form with wide and flat bottom to avoid tipping

Available in 15 different shapes and sizes

FDA food grade compliant material

Bisphenol A (BPA) and Latex-Free





DISPOSABLE WEIGH BOATS & CANOES

PICTURE	DESCRIPTION	SIZE	COLOR	QUANTITY	PRODUCT CODE	DIMENSIONS	CAPACITY
 	Weigh Boat, Antistatic	Small	Transparent Light White	100	6032-01	1.6 x 1.6 x 0.3" (41 x 41 x 8 mm)	0.25 oz. (7 mL)
				500	6032-02		
				2500	6032-03		
			Black	100	6041-01	1.6 x 1.6 x 0.3" (41 x 41 x 8 mm)	0.25 oz. (7 mL)
				500	6041-02		
				2500	6041-03		
		Medium	Transparent Light White	100	6033-01	3.5 x 3.5 x 1" (89 x 89 x 25 mm)	3.52 oz. (100 mL)
				500	6033-02		
				2500	6033-03		
			Black	100	6042-01	3.5 x 3.5 x 1" (89 x 89 x 25 mm)	3.52 oz. (100 mL)
				500	6042-02		
				2500	6042-03		
		Large	Transparent Light White	100	6034-01	5.5 x 5.5 x 1" (140 x 140 x 25 mm)	11.61 oz. (330 mL)
				500	6034-02		
				2500	6034-03		
Black	100		6043-01	5.5 x 5.5 x 1" (140 x 140 x 25 mm)	11.61 oz. (330 mL)		
	500		6043-02				
	2500		6043-03				
	Weigh Boat, Diamond, Antistatic	Medium	Transparent Light White	250	6039-01	3.1 x 2.4 x 0.7" (80 x 60 x 18 mm)	1.06 oz. (30 mL)
				500	6039-02		
	Weigh Boat, Hexagonal, Antistatic	Medium	Transparent Light White	250	6040-01	OD: 3", Dep: 3/4" (OD: 76 mm, Dep: 19 mm)	1.76 oz. (50 mL)
				500	6040-02		
 	Weigh Canoe, Antistatic	Small	Transparent Light White	100	6035-01	2.1 x 1.5 x 0.5" (53 x 38 x 13 mm)	0.53 oz. (15 mL)
				500	6035-02		
				100	6044-01		
			Black	500	6044-02	2.1 x 1.5 x 0.5" (53 x 38 x 13 mm)	0.53 oz. (15 mL)
				250	6036-01		
				500	6036-02		
		Medium	Transparent Light White	250	6036-01	2.7 x 2 x 0.6" (68 x 50 x 16 mm)	0.88 oz. (25 mL)
				500	6036-02		
			Black	250	6045-01	2.7 x 2 x 0.6" (68 x 50 x 16 mm)	0.88 oz. (25 mL)
				500	6045-02		
		Large	Transparent Light White	100	6037-01	5.2 x 3.2 x 1" (133 x 85 x 26 mm)	4.22 oz. (120 mL)
				500	6037-02		
100	6046-01						
Black	500		6046-01	5.2 x 3.2 x 1" (133 x 85 x 26 mm)	4.22 oz. (120 mL)		
	100		6046-01				
	500		6046-02				
X-Large	Transparent Light White	250	6038-01	7.1 x 4.3 x 1" (181 x 110 x 26 mm)	8.45 oz. (240 mL)		
		500	6038-02				

FOR MORE INFORMATION

USA 1-800-932-1039 | CAN 1-800-665-6334 | AUS 1300 786 392 | WWW.MEDISCA.COM



DISCLAIMER: The information presented is not geared toward any indication whatsoever and is solely at the discretion and judgment of the duly licensed medical practitioner. It is the responsibility of the licensed pharmacist or other appropriately licensed professional to determine if the products referenced herein are appropriate for the specific use. MEDISCA cannot be responsible and shall not be responsible for customer's use of the information contained in this suggested product information sheet.