



PRODUCT NO.: 6041







Weigh Boat, Anti-Static, 1.6 x 1.6 x 0.3", Small, Black

Description: Anti-static disposable polystyrene weight boats are flexible for handling liquids and solids. Square with stable rounded bottom, smooth surface and contoured sides.

*Weigh boats contain carbon black pigment and are not recommended for use with oils.

FEATURES:

-  Easy pour spout
-  Boats are shallow with a wide, flat bottom to resist tipping
-  Rounded corners simplify transfer
-  Product is resistant to dilute and weak acids, aqueous solutions, alcohols, and bases

TECHNICAL SPECIFICATIONS:

Capacity (mL)	7
Dimension (mm/in)	41 x 41 x 8/1.6 x 1.6 x 0.3
Color	Black
Material	Polystyrene (PS), Anti-Static

USA | CANADA | AUSTRALIA PACIFIC | EUROPE
WWW.MEDISCA.COM