



**PRODUCT NO.: 6043**







## **Weigh Boat, Anti-Static, 5.5 x 5.5 x 1", Large, Black**

**Description:** Anti-static disposable polystyrene weight boats are flexible for handling liquids and solids. Square with stable rounded bottom, smooth surface and contoured sides.

\*Weigh boats contain carbon black pigment and are not recommended for use with oils.

### **FEATURES:**

-  Easy pour spout
-  Boats are shallow with a wide, flat bottom to resist tipping
-  Rounded corners simplify transfer
-  Product is resistant to dilute and weak acids, aqueous solutions, alcohols, and bases

### **TECHNICAL SPECIFICATIONS:**

Capacity (mL)	330
Dimension (mm/in)	140 x 140 x 25/5.5 x 5.5 x 1
Color	Black
Material	Polystyrene (PS), Anti-Static

USA | CANADA | AUSTRALIA PACIFIC | EUROPE  
WWW.MEDISCA.COM