PRODUCT NO.: 6046





## Weigh Canoe, Anti-Static, 5.2 x 3.2 x 1", Large, Black

**DESCRIPTION:** This weigh canoe protects the pans of weighing balances from damage and provides a convenient way to transfer weighed materials to another container.

## FEATURES:

- 2 Pour-spout design easily dispenses samples without flexing the dish
- Z Flat-bottom dishes are from 12-mil thick polystyrene for better strength, stability and handling

## **TECHNICAL SPECIFICATIONS:**

Capacity (mL)	120
Dimension (mm/in)	133 x 85 x 26/5.2 x 3.2 x 1
Color	Black
Material	Polystyrene (PS), Anti-Static