



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6058





Hood Accessories, 3' Base Stand

w/Casters 36" Working Height (for SafeAiRx)

DESCRIPTION: Constructed of 14-gauge epoxy coated steel, these custom carts are designed exclusively for our ductless hoods and workstations. All carts have a height of 36" and feature a two-position movable base bar for wheelchair or lab chair compatibility.

APPLICABLE PRODUCTS:

-  6583: 3' Ductless Fume Hood, Particulate & Vapor Extraction Enclosure (SafeAiRx)
-  6059: Hood, 3' SafeAirX Double HEPA Hood w/Canopy & Stainless Steel Work Surface, 115V/60Hz

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM