




**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6061**



## Tray Holder Insert (Labconco)




Type 304 Stainless Steel

**DESCRIPTION:** Pins hold a variety of wider-mouth flasks, beakers and other labware in inverted position. Type 304 stainless steel pins have plastic tips to protect glassware from scratching. Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

### TECHNICAL SPECIFICATIONS:

Dimensions, W x D x H (in)	19.3 x 8.8 x 8.0
Estimated Shipping Weight (kg/lbs)	1.1/2.5
Rack Location	Lower, Upper

### APPLICABLE PRODUCTS:

-  6066-01: SteamScrubber Glassware Washer, Labconco (W 24.2" x D 27.5" x H 34.2" – 36.2", Freestanding)
-  7496-01: Steamscrubber Glassware Washer, Labconco (W 24.2" x D 27.5" x H 34.2" – 36.2", Freestanding)
-  8654-01: Steamscrubber Glassware Washer, Labconco (W 24.1" x D 27.4" x H 34.1" – 36.1", Undercounter)

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)