



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6063



10-Pin Labware Insert (Labconco)

Type 304 Stainless Steel

DESCRIPTION: Pins hold a variety of wider-mouth flasks, beakers and other labware in inverted position. Type 304 stainless steel pins have plastic tips to protect glassware from scratching. Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

TECHNICAL SPECIFICATIONS:

Dimensions, W x D x H (in)	19.3 x 8.8 x 8.0
Estimated Shipping Weight (kg/lbs)	1.1/2.5
Rack Location	Lower, Upper

APPLICABLE PRODUCTS:

- 6066-01: SteamScrubber Glassware Washer, Labconco (W 24.2" x D 27.5" x H 34.2" – 36.2", Freestanding)
- 7496-01: Steamscrubber Glassware Washer, Labconco (W 24.2" x D 27.5" x H 34.2" – 36.2", Freestanding)
- 8654-01: Steamscrubber Glassware Washer, Labconco (W 24.1" x D 27.4" x H 34.1" – 36.1", Undercounter)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM