



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6064**



## **Work Surface Plate, Labconco** (4 ft, Deep, Black Epoxy)

**DESCRIPTION:** Black chemical-resistant work surface is contoured to fit the dimensions of the enclosure and to contain spills. Can be used for XPert/RXPert Filtered Balance System & Purifier Vertical Clean Bench.

### **TECHNICAL SPECIFICATIONS:**

Dimensions, W x D x H (in/cm)	48 x 35.5 x 1/122 x 90 x 2.5
Shipping Weight (lbs/kg)	125/56

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)