



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6111



Mixer, Mazerustar KK-400W, Counter Balance Weight Set







Description: This new counter balance weight set can be placed directly in the Mazerustar KK400W holder facilitating the counter balance step required before initiating any cycle on the Maz. The corresponding mass engraved of each weight facilitates identification

TECHNICAL SPECIFICATIONS:

Total Weight (g)	430
Material	Acetal
Color	White

* Please note there is a margin of error of $\pm 12g$

INCLUDES:

-  1 x 200g
-  1 x 100g
-  1 x 50g
-  1 x 40g
-  1 x 25g
-  1 x 15g

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM