




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6117



Mixer, Mazerustar KK-400W, Counter Balance Weight, 15g

Description: These new counter balance weights can be placed directly in the Mazerustar KK400W holder facilitating the counter balance step required before initiating any cycle on the Maz. The corresponding mass engraved of each weight facilitates identification.

TECHNICAL SPECIFICATIONS:

Total Weight (g)	15
Material	Acetal
Color	White

* Please note there is a margin of error of $\pm 2g$

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM