



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING



**PRODUCT NO.: 6129**



## Mixer, Mazerustar KK-400W Spacer for Mixing Adapters

**DESCRIPTION:** This spacer is designed to help support the adapter in the Maz. This spacer is placed below the mixing adapter being used in the Mazerustar KK-400W.

### FEATURES:

-  Extends the life of your adapters
-  Compatible with MD adapters

### TECHNICAL SPECIFICATIONS:

Height (mm)	27
Diameter (mm)	77
Weight (g)	118

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)