



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING




PRODUCT NO.: 6145



Stopper, 20 mm, Snap-On, Chlorobutyl, Grey (Sterile)

DESCRIPTION: The snap-on stoppers lock into the vial, providing a tight seal and keeping the stopper from backing out.

FEATURES:

-  Pliable elastomer material
-  No fray, no tear design
-  Air tight seal

TECHNICAL SPECIFICATIONS:

Inner Diameter (mm)	13.12
Outer Diameter (mm)	18.80
Thickness (mm)	3.28

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM