



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6165**



**BOTTLE ADAPTER (32 – 33 mm)**  
ETFE, For BrandTech Bottletop Dispensers

**DESCRIPTION:** This bottle adapter is made of Ethylene TetraFluoroEthylene (ETFE). It should be used in place of the standard Polypropylene (PP) adapters that come with the Bottletop dispensers, in cases where higher chemical resistance is required.

**TECHNICAL SPECIFICATIONS:**

Outer Thread Dimension, For Dispenser (mm)	45
Inner Thread Dimension, For Bottle (mm)	32 – 33
Material	ETFE