




Vials, Open, Amber Glass, 100 mL, 20 mm (Sterile)

Description: All vials are sterile. Lot release sterility testing conforms to requirements set forth in USP <71>. Stopper and seal are not included.

TECHNICAL SPECIFICATIONS:

Filling Volume (mL)	106.5 ± 2.0
Height (mm)	94.5 ± 0.5
Thickness of bottle-lip (mm)	4.2 ± 0.2
Height of bottle-neck (mm)	3.5 ± 0.3
Effective height of bottle-body (mm)	62.0 ± 2.0
Outer diameter of bottle-mouth (mm)	20.0 ± 0.2
Inner diameter of bottle-mouth (mm)	12.6 ± 0.2
Outer diameter of bottle-neck (mm)	16.7 ± 0.3
Body diameter (mm)	51.6 ± 0.3
Material	Borosilicate glass



APPLICABLE PRODUCTS:

- 🚫 6118: Seals, Aluminum, Center Tear-Out, 13 mm, Silver (Sterile)
- 🚫 6119: Stopper, 20 mm, 2 Prongs, Chlorobutyl, Grey (Sterile)
- 🚫 6141: Seals, Aluminum, Flip-Off, 20 mm, Blank Top, Red (Sterile)
- 🚫 6142: Seals, Aluminum, Flip-Off, 20 mm, Blank Top, Blue (Sterile)
- 🚫 6143: Seals, Aluminum, Flip-Off, 20 mm, Blank Top Green (Sterile)
- 🚫 6145: Stopper, 20 mm, Snap-On, Chlorobutyl, Grey (Sterile)
- 🚫 6146: Stopper, 13 mm, Snap-On, Chlorobutyl, Grey (Sterile)

USA | Canada | Australia | International
medisca.com

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.