




Vials, Open, Amber Glass, 30 mL, 20 mm (Sterile)

Description: All vials are sterile. Lot release sterility testing conforms to requirements set forth in USP <71>. Stopper and seal are not included.

TECHNICAL SPECIFICATIONS:

| | |
|--------------------------------------|--------------------|
| Filling Volume (mL) | 33.0 ± 1.0 |
| Height (mm) | 62.8 ± 0.5 |
| Thickness of bottle-lip (mm) | 4.2 ± 0.2 |
| Height of bottle-neck (mm) | 3.5 ± 0.3 |
| Effective height of bottle-body (mm) | 40.0 ± 1.0 |
| Outer diameter of bottle-mouth (mm) | 19.8 ± 0.2 |
| Inner diameter of bottle-mouth (mm) | 12.6 ± 0.2 |
| Outer diameter of bottle-neck (mm) | 16.7 ± 0.3 |
| Body diameter (mm) | 35.8 ± 0.3 |
| Material | Borosilicate glass |



APPLICABLE PRODUCTS:

- 📄 6118: Seals, Aluminum, Center Tear-Out, 13 mm, Silver (Sterile)
- 📄 6119: Stopper, 20 mm, 2 Prongs, Chlorobutyl, Grey (Sterile)
- 📄 6141: Seals, Aluminum, Flip-Off, 20 mm, Blank Top, Red (Sterile)
- 📄 6142: Seals, Aluminum, Flip-Off, 20 mm, Blank Top, Blue (Sterile)
- 📄 6143: Seals, Aluminum, Flip-Off, 20 mm, Blank Top Green (Sterile)
- 📄 6145: Stopper, 20 mm, Snap-On, Chlorobutyl, Grey (Sterile)
- 📄 6146: Stopper, 13 mm, Snap-On, Chlorobutyl, Grey (Sterile)

USA | Canada | Australia | International
medisca.com

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.