


PRODUCT NO.: 6205




Vials, Open, Amber Glass, 2 mL, 13 mm (Sterile)

Description: All vials are sterile. Lot release sterility testing conforms to requirements set forth in USP <71>. Stopper and seal are not included.

TECHNICAL SPECIFICATIONS:

Filling Volume (mL)	3.5 ± 0.5
Height (mm)	35.0 ± 0.5
Thickness of bottle-lip (mm)	3.9 ± 0.2
Height of bottle-neck (mm)	4.0 ± 0.3
Effective height of bottle-body (mm)	23.0 ± 1.0
Outer diameter of bottle-mouth (mm)	13.0 ± 0.2
Inner diameter of bottle-mouth (mm)	7.0 ± 0.2
Outer diameter of bottle-neck (mm)	10.0 ± 0.3
Body diameter (mm)	16.0 ± 0.15
Material	Borosilicate glass



APPLICABLE PRODUCTS:

- 6118: Seals, Aluminum, Center Tear-Out, 13 mm, Silver (Sterile)
- 6119: Stopper, 20 mm, 2 Prongs, Chlorobutyl, Grey (Sterile)
- 6141: Seals, Aluminum, Flip-Off, 20 mm, Blank Top, Red (Sterile)
- 6142: Seals, Aluminum, Flip-Off, 20 mm, Blank Top, Blue (Sterile)
- 6143: Seals, Aluminum, Flip-Off, 20 mm, Blank Top Green (Sterile)
- 6145: Stopper, 20 mm, Snap-On, Chlorobutyl, Grey (Sterile)
- 6146: Stopper, 13 mm, Snap-On, Chlorobutyl, Grey (Sterile)

USA | Canada | Australia | International

medisca.com