




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6217



LIFTING HANDLE EXTENSION

For Profiller 3700 Capsule Machine

DESCRIPTION: This extension helps to move the lifting handle, of the Profiller 3700 capsule machine, closer to the operator. This makes it easy to lock the capsules when used under a hood. It also modifies the reach distance to separate or lock the capsules reducing the level of discomfort to the operator.

Included in Kit: LHE Bracket (1), Bolt (1), Nut (1)

Required: A 17 mm wrench is required for installation