



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6226**



Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.

## **MOBILE CASTER KIT, LABCONCO**

For Hydraulic Lift Base Stands

**DESCRIPTION:** This caster kit adds mobility to the hydraulic lift base stands. It replaces the fixed feet that come with the base stand and is quiet in operation.

**Included in Kit:** 4 non-marking Polyurethane Casters with Tie Locks (4)

### **TECHNICAL SPECIFICATIONS:**

Dimensions, Ø (cm/ in)	12.7 / 5
Shipping Weight (Kg / lb)	5 / 10
Load Rating per Caster (lb)	250