

PRODUCT NO.: 6242



Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.

Chemotherapy Spill Kit

Description: This chemotherapy spill kit is designed for the cleanup of spills involving cytotoxic drugs. It helps to solidify and contain hazardous chemical spills for quick and easy cleanup and disposal with minimal exposure.

Included in Package: Nitrile Gloves (2 Pairs), Nitrile Chemotherapy Gloves (2 Pairs), Full-Face Safety Shield (1), Chemotherapy Gown (1), Green Z Solidifier (1x 21 g), Scoop (1), Scraper (1), Wiper Pad (2), Non-slip Shoe Cover (1), Red Bio-hazard Waste Bag (1), Twist Tie (1), Identification Tag (1), Instruction (1)

Additional Information : Green Z Solidifier 1oz solidifies approximately 1L.

USA | Canada | Australia | International

medisca.com