

PRODUCT NO.: 6256



Bottle, Natural, Dropper 1 oz (30 mL), w/Cap

Description: This bottle is squeezable and can be used for nasal applications. It can also be used as a dropper bottle.

Technical Specifications:

| | | |
|------------------------------|---------|--------------------------|
| Capacity (oz / mL) | | 1 / ~30 |
| Bottle Dimension, Ø × H (in) | | 1.236 × 2.950 ± 0.05 |
| Neck Finish | | 18-410 |
| Overflow capacity (oz / mL) | | 1.27 ± 0.08 / 37.7 ± 2.5 |
| Color | | Natural |
| Bottle Type | | Boston Round |
| Material | Cap | PP |
| | Dropper | LDPE |
| | Bottle | LDPE |

***Note:** Not intended as medical device

USA | Canada | Australia | International
medisca.com

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.