



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6283**





## **Hood Accessories, 4' Base Stand**

w/Casters 36" Working Height (for SafeAirRx)

**DESCRIPTION:** Constructed of 14-gauge epoxy coated steel, these custom carts are designed exclusively for our ductless hoods and workstations. All carts have a height of 36" and feature a two-position movable base bar for wheelchair or lab chair compatibility.

### **APPLICABLE PRODUCTS:**

-  6571: 4' Ductless Fume Hood, Particulate & Vapor Extraction Enclosure (SafeAirRx)
-  6057: Hood, 4' SafeAirX Double HEPA Hood w/Canopy & Stainless Steel Work Surface, 115V/60Hz

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)