



Capsule Machine Accessories, Profiller 3700 Sliding Sheet Set, Size 3

DESCRIPTION: Replacement body sheet set for the Profiller 3700, Size 3. This is required when filling capsules for size 3. This is typically included in a combo or change part for this corresponding size.

USA | CANADA | AUSTRALIA PACIFIC | EUROPE WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.