



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6300



Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.

## **DRAFT SHIELD, OHAUS**

### **For Adventurer Analytical Balance**

**DESCRIPTION:** This draft shield helps to prevent unstable readings due to air currents in the laboratory or inside the hood. It also reduces the footprint, of the balance, when fully opened.

The two-door configuration, allows the draft shield to open wide, to give enough space to load samples during weighing procedures.

#### **TECHNICAL SPECIFICATION:**

Material of Construction	Glass with Plastic Lining and handles
--------------------------	---------------------------------------